



96x DIMPLES

0.25

R0.125

**DETAIL A**

SCALE 2 : 1

UNLESS OTHERWISE SPECIFIED: DIMENSIONS ARE IN INCHES TOLERANCES: ONE DECIMAL: 0.05 TWO PLACE DECIMAL: 0.01 THREE PLACE DECIMAL: 0.005		PROPRIETARY AND CONFIDENTIAL: THE INFORMATION CONTAINED IN THIS DRAWING IS THE SOLE PROPERTY OF V&P SCIENTIFIC, INC. ANY REPRODUCTION IN PART OR AS A WHOLE WITHOUT THE WRITTEN PERMISSION OF V&P SCIENTIFIC, INC. IS PROHIBITED.		DEBURR AND BREAK SHARP EDGES		DO NOT SCALE DRAWING		REVISION	
						<b>V&amp;P SCIENTIFIC, INC.</b>			
						TITLE: <b>VP 530D-PTFE-1</b>			
						MATERIAL: <b>PTFE</b>		DWG NO. <b>498</b>	
						FINISH:		SCALE:1:1	
								SHEET 1 OF 1	

	NAME	SIGNATURE	DATE		
DRAWN	AL		11/11/21		
CHK'D					
APPV'D					
MFG					
Q.A					